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WORKING PROCEDURE FOR CHROMATOGRAPHIC COLUMN (introduction of the sample by dry method):

1. For a chromatographic column 42 cm high and 2.5 cm in diameter (reaction crude on the order of 500 mg), we take about 70 cm³ of silica-60 for chromatography measured in a measuring cylinder.

For a chromatographic column 30 cm high and 1.5 cm in diameter (reaction crude on the order of 250 mg), we take about 25 cm³ of silica-60 for chromatography measured in a measuring cylinder.

2. Prepare the gel with silica-60 and the solvent or mixture of solvents that we will use as the eluent.

3. Fill the column with the eluent up to a fifth or a quarter of its height.

4. Put the gel into the column.

5. Allow some time for the gel to settle. Then, we speed up the process by opening the key a little. The eluent must be allowed to reach the top of where the gel reaches.

6. The mixture to be separated is dissolved in the minimum volume of the major solvent to be used in the eluent in a 100 mL Erlenmeyer flask with ground glass joint. For the large chromatographic column (42x2.5cm) 2.5 cm³ of silica-40 is added and for the small chromatographic column (30x1.5cm) 1.2 cm³ of silica-40 is added.

7. The suspension is concentrated to dryness in the rotavapor.

8. Silica-40 with the absorbed sample is poured cautiously with the help of a spoon-shaped spatula and a funnel for solids inside the column.

9. On top of the above solid, a layer of about 2 millimeters of clean sea sand is added with a spoon-shaped spatula and a funnel for solids .

10. Add the eluent carefully and through the wall of the column and at the same time open the key to introduce the sample into the stationary phase.

11. Add the eluent until filling the column.

The level of the eluent should not be allowed to be below the upper level of the gel.

It is convenient to filter the fractions that are collected, since part of the stationary phase is not retained by the filter plate of the column.

Work in the fume chamber. Chlorinated solvents are harmful, and when handling silica-60, particles remain suspended in the air.

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